





**RAMPAS E PATAMARES** 

|               | 1      | 2      | 3      | 4      | 5      |
|---------------|--------|--------|--------|--------|--------|
| Set Point     | 150 °C | 200 °C | 250 °C | 300 °C | 350 °C |
| Tempo Rampa   | 15 MIN | 20 MIN | 25 MIN | 30 MIN | 35 MIN |
| Tempo Patamar | 30 MIN | 35 MIN | 40 MIN | 45 MIN | 50 MIN |

|   |   |
|---|---|
| Patamar                                   | Numero Ciclos                           |
| <input checked="" type="checkbox"/> Tempo | 10                                      |
| <input type="checkbox"/> Temperatura      | Fin de Ciclo:                           |
| Faixa Sup: 20 °C                          | <input type="checkbox"/> Controle       |
| Faixa Inf: 30 °C                          | <input checked="" type="checkbox"/> OFF |

**RAMPAS E PATAMARES** 

☒ SEGUIMENTO 1

|               |        |
|---------------|--------|
| Set Point     | 150 °C |
| Tempo Rampa   | 15 MIN |
| Tempo Patamar | 30 MIN |

# ***DIGESTOR BLOCK WITH TOUCH SCREEN RAMPS AND LANDERS***

## ***TE-041/25***

Used to digest the most diverse types of samples, such as plants, food, rubber, among others, for later analysis of nitrogen / protein.

## Technical Characteristics

### TE-041/25

- Temperature: Ambient +7°C to 450°C;
- Temperature Controller: FE15J Digital controller board with PID system and ramps and steps with 1 program divided into 5 segments (temperatures) selectable with a maximum of 5 ramps and 5 steps (if 5 temperatures are selected). Possibility of 10 cycle repetitions;
- Selection of counting type: By time (counts the time regardless of temperature and obeys the programmed time) or by temperature (only starts counting if the system is within the programmed limit);
- Control selection after ramps and soaks: Turns off the control or controls at the last set point;
- Operation screen: TOUCH SCREEN 4.3;
- Sensor: J type with stainless spring;
- Accuracy:  $\pm 1^{\circ}\text{C}$ ;
- Uniformity:  $\pm 3^{\circ}\text{C}$ ;
- Safety: Shielded resistance avoiding contact with sulfuric acid and type B circuit breaker;
- Block: In die-cast aluminum with 45 mm hole depth;
- Cabinet: Stainless steel 304;
- Dimensions: W=320 x D=430 x H=330 mm;
- Weight: 18 kg;
- Power: 2200 Watts;
- Voltage: 220 Volts;
- ACCOMPANIES: - 01 gallery in 304 stainless steel - 01 digital temperature controller with ramps and soaks - 40 micro tubes in 90 ml borosilicate glass  $\varnothing 25 \times 250$  mm - Instruction Manual with Warranty Term;

## Benefits and Advantages

- Compact equipment
- User-friendly touch screen display
- It has date and time on the display
- It has Stand-by mode
- After starting the process, it takes place automatically, providing agility
- It can have 5 segments, therefore 5 ramps and 5 steps
- Allows control of the time that the sample must remain at a certain temperature (level)
- Allows control of acceptable temperature variation in the process
- The ramp, when well dimensioned, allows the temperature rise to be smooth, resulting in low overshoot and better temperature distribution across the block
- At the end of the fulfillment of the programmed segments there is a resource item called looping (cycles), it can be 0 to 10, providing agility
- Allows the jumper option, in which two thresholds with the same set point can be programmed, providing plateau with longer time
- At the end of the process there is the control option (remaining at constant temperature) or the end option, the client choice
- It has a controller box separate from the block, which can be left outside the hood, resulting in a longer useful life to the equipment for being more protected from acidic gases and vapors
- Possibility of adaptations according to the client needs, makes the equipment already of line a special equipment.