

+55 (19) 2105-6161 comex@tecnal.com.br





## LABORATORY Y-BLENDER TE-201/10

Used to homogenize samples into powder or granules.



João Leonardo Fustaino, nº 325 Distrito Industrial Uninorte Piracicaba/SP-Brasil • CEP 13.413-102



## **Technical Characteristics**

## TE-201/10

- Shaking Amplitude: 360°;
- Homogenizer: constructed of 304 stainless steel with polished finish, internal and external;
- Motor: 1/3 three-phase with reducer;
- Covers: polished 304 stainless steel, removable, with silicone seal:
- Total volume: 18 Liters;
- Gabinet: in carbon steel with anti-corrosion treatment and electrostatic painting;
- Overall dimension: L = 830 X D=420 x H=750 mm; Usable volume: 7,5 liters;
- Weight: 33 kg;
- Power: 250 Watts;

- Voltage: 220V +/-5% 50/60 Hz;
- Include: 02 Extra fuses Instruction Manual and warranty term;
- NOTE: The informed volume has as reference product with density 0,7Kg/dm3;
- Rotation range: 10 to 40 RPM ;
- Accuracy: ±2 RPM ;
- Rotation control: Analog via Potentiometer and frequency variator;
- RPM reading: Via digital tachometer;

## **Benefits and Advantages**

- · Easy-to-remove cap at three ends for cleaning, sample input and material output
- Vat dyes in stainless steel 304 with welding and fine polishing to have a good asepsis avoiding contamination and for longer service life
- Brushless induction Motor, can be used for long periods
- Strict quality control, in which checks and tests ensure the perfect operation of the equipment, providing safety and customer satisfaction
- Customer service, to answer questions and provide explanations about the equipment and methodologies
- Possibility of adaptations according to the needs of the customer, makes the line equipment a special equipment

