



# OPEN-CELL STIRRER

## TE-083

Used to analyze disintegrated sugar cane cells.

## Technical Characteristics

### TE-083

- Rotation: Fixed at 60 RPM;
- Accuracy:  $\pm 1$  RPM;
- Agitation range: 360°;
- Motor: Reducer 84 Watts;
- Capacity: 02 containers of 4 liters;
- Fixing of the containers: Through handles;
- Cabinet: In carbon steel with anti-corrosive treatment and electrostatic painting;
- Dimensions: W = 520 x D = 500 x H = 700 mm;
- Weight: 20 kg;
- Power: 100 Watts;
- Voltage: 220 Volts;
- Include: - 02 Plastic containers (4 liters) - 02 extra fuses - Instruction Manual with Warranty Term;
- Optional: Protection cabinet with emergency stop sensor;

### Benefits and Advantages